



COLPOSCOPES

Take a look at our new C-100 series and discover why it is your best option

COLPOSCOPES



ECLERIS 
INNOVATION VALUE



info@ecleris.com
www.ecleris.com

COLPOSCOPE

C-100 SERIES

- **HIGH IMAGE QUALITY**

High quality wide field optics, 5 levels of magnification. Integrated fine focus objective lens and fully adjustable dual maneuvering handles.

- **OUTSTANDING MOBILITY**

All rotary components are assembled on bearings providing smooth and stable motion.

- Optional monitor Stand

- Optional cameras :
Procam standard and
high definition resolution

- New led 80w light source
- Universal connection for all brands of fiber optic cables
- Electronic dimmer
- 50.000 hours life LED

- **Reliable support system**

Telescopic compact option

Mobile base with rubber wheels for firm and stable movement with fixing position breaks

