





## “Talking” Automatic Blood Pressure Monitor

Model: CF175 



- Talking Function  
multi language (EN/RU/FR/AR/FA/TR)
- Cuff Detection
- Movement Detection
- Real Fuzzy Technology
- Hypertension Risk Indication
- Irregular Heartbeat Detection (IHB)
- 90 memories
- Average of Last 3 Readings
- One-touch Automatic Operation
- Latex-Free Patented Universal Cone Cuff
- Adaptor Socket
- Lifetime Calibration
- Carrying bag



### Clinically Accurate

This blood pressure monitor has passed the criteria of the European Society of Hypertension by achieving all requirements. This monitor further meets the AAMI accuracy standard requirements (AAMI/ANSI/ISO 81060-2)



Universal Cone cuff fits arm circumference:  
24~40cm (9.4"~15.7")  
Other size cuffs available upon request

Model	Qty per carton	Carton volume
CF175f	18 pcs	0.068 cbm/ctn
QTY per 20'	QTY per 40'	QTY per 40'HQ
7506 pcs	15030 pcs	17460 pcs

● Storage and Transportation Condition  
-10°C~60°C(14°F~140°F)  
10%~90% RH, 700~1060 hPa

Patent Nos:  
US: 7597667B2; CN: 469964  
TW: 191076; TW: 1236891; TW: 1322262  
Other patents pending

DISTRIBUTOR:  
**Beracah**  
MEDIC A.MX  
ventas@beracahmedica.mx

### Real Fuzzy Technology

Real Fuzzy Technology determines ideal cuff pressure based on one's systolic blood pressure and arm size. Users no longer need to pre-set the inflation level before measurement. The technology eliminates arm discomfort through unnecessary high inflation settings and also prevents inaccurate readings caused by erroneous cuff inflation levels.