

# Vortex-Genie® Mixers

THE ORIGINAL **Beracah**  
MEDICA



## Vortex-Genie® 1

- Violent single high-speed vortexing
- Easy-to-use Touch On operation
- Large, 6mm orbit for aggressive vortexing of difficult samples
- May be used with Pop-off Cup (shown) or optional 3-inch Platform
- May be used in cold rooms and incubators

## Vortex-Genie® 2

- OUR MOST POPULAR MODEL
- Variable speed for precise mixing from gentle to vigorous
- Broad range of attachments for most tubes, plates and other containers
- Hands-Free or Touch On control
- May be used in cold rooms and incubators

### SPECIFICATIONS

Vortex-Genie Model	Vortex-Genie 1	Vortex-Genie 2
Supplied	Pop-off Cup	3-inch Platform, Pop-off Cup
Speed (RPM)	3200 (2700 for 50Hz models)	600-3200 (2700 for 50Hz models)
Timer - Touch On: Hands-Free:	N/A	N/A
Dimensions (base) (D x W x H)	165 x 122 x 165mm (6.5 x 4.8 x 6.5in)	165 x 122 x 165mm (6.5 x 4.8 x 6.5in)
Weight	4 Kg (8.8 lb)	4 Kg (8.8 lb)
Markings / Warranty	ETL, CE / 2-Years	ETL, CE / 2-Years

### CATALOG NUMBERS

120V	SI-0136	SI-0236 (Model G560)
230V - No Plug	SI-0146	SI-0246 (Model G560E)
230V - European Plug	SI-0156	SI-0256 (Model G560E)
230V - British Plug	SI-0166	SI-0266 (Model G560E)
230V - Swiss Plug	SI-0176	SI-0276 (Model G560E)
240V - Australian Plug	SI-0196	SI-0297 (Model G560E)
100V	SI-0186	SI-0286

See page 3 for accessories.

U.S. Pats. 4,781,487  
5,707,861